# Dr. Huang Lab Biological Safety Task Procedure Human/Mouse Material

Date: May 1, 2013

Author: Zhaojun Yin Phone: 517-897-3008

Job Task

Human/mouse Cell Line Culture

# **Description**

Any procedure involved in the manipulation of human/mouse cell lines for culturing in this laboratory. May include but not limited to: Transferring, splitting, feeding, harvesting, inoculation/treating, or freezing of the cells. Opening cell culture flasks for any manipulations. Adding or removing material from flasks, dishes and vials while using aseptic technique.

#### Hazards

Bloodborne Pathogens, Bleach or other Disinfectant

# **Engineering Controls**

Biological Safety Cabinet Eyewash Sink Mechanical Pipettor

# **Personal Protective Equipment**

Gloves Lab Coat

### **Work Practices**

- ✓ Turn on cabinet at least 10 minutes prior to use.
- ✓ Before beginning work, close culture room door.
- ✓ Spray all items down with at least 70% ETOH or other disinfectant prior to bringing into cabinet.
- ✓ Add disinfectant (bleach) to vacuum flask prior to use.
- ✓ Open items in biosafety cabinet only.
- ✓ Use aseptic technique and work from clean to dirty inside cabinet.
- ✓ Solid waste must be in a biohazard bag (2 with serological pipettes), a 2nd container that is leak-proof, puncture resistant and closed when unattended.
- ✓ Dispose of biohazardous waste daily: liquid (vacuum flasks) or solid.
- ✓ Liquid waste must be either treated by chemical decontamination or autoclaved immediately after collection and prior to disposing by pouring into sanitary sewer.
- ✓ Solid waste must be autoclaved and then placed in opaque bag into dumpster for disposal.
- ✓ Spray all items down with 10% Bleach, 70% ETOH, or other disinfectant prior to removing from cabinet.
- ✓ Clean cabinet by removing any waste or other items. Do not store items in cabinet or on grills.
- ✓ Disinfect cabinet with 10% bleach or (other appropriate disinfectant like Rocal) as directed by label. Rinse disinfectant if needed by using 70% ETOH or water.
- ✓ Remove gloves, dispose of them in solid biohazard waste, and wash hands with soap and water.